

Anti-Human BACE (C) Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01254

3th Edition

Description: The current understanding of the role of amyloid in Alzheimer's disease (AD) has been established by a remarkable congruence of multiple disciplines: neuropathology, biochemistry, molecular biology, and epidemiological genetics.

Antigen: Synthetic peptide of C terminal part of human BACE Protein (CLRQQHDDFADDISLLK)

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Cross reacts with both mouse and rat.