

**IgG (H+L) Anti-Rabbit Goat IgG Fab'- HRP**

**Polyclonal Antibody**

**Cat.NO.: PA01259**

---

3th Edition

**Description:** IgG (H+L) Anti-Rabbit Goat IgG Fab'- HRP, 0.2 mg. For research use only, not for use in diagnostic procedures.

**Antigen:** Rabbit IgG (H+L)

**Form:** Lyophilized product from 1% BSA in PB (pH6.5) (containing ProClin®300 as preservative)

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** The antibody is absorbed by goat IgG, mouse IgG, human IgG and rat IgG.