

IgG (H+L) Anti-Rabbit Goat IgG Fab'- HRP

Polyclonal Antibody

Cat.NO.: PA01259

3th Edition

Description: IgG (H+L) Anti-Rabbit Goat IgG Fab'- HRP, 0.2 mg. For research use only, not for use in diagnostic procedures.

Antigen: Rabbit IgG (H+L)

Form: Lyophilized product from 1% BSA in PB (pH6.5) (containing ProClin®300 as preservative)

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: The antibody is absorbed by goat IgG, mouse IgG, human IgG and rat IgG.