

Kappa Light Chain Goat IgG Biotin (Monkey Absorbed)

Polyclonal Antibody

Cat.NO.: PA01273

3th Edition

Description: Anti-Human IgG (goat IgG) is often used as a secondary antibody for the detection of Human IgG. Following recent advances in the development of antibody-based therapy, such as the development of Human antibodies, this secondary antibody has been used for the purpose of investigating metabolism in the living body following administration of this antibody.

Antigen: Purified Human Kappa light chain

Form: Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: