

LRG (329) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01293

3th Edition

Description: Among the diseases which should be distinguished from iNPH (idiopathic normal pressure hydrocephalus) of elderly people showing symptoms such as gait disorder and dementia, there are some neurodegenerative diseases such as Alzheimer's disease, FTLD (frontotemporal lobar degeneration) and disorder associated with Parkinson's disease.

Antigen: Synthetic peptide of the C terminal part of Human LRG (CAGPEAVKQGQTLLAVAKSQ)

Form: Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: