

**Prion Protein (N) Anti-Human Rabbit IgG Affinity Purify
Polyclonal Antibody**

Cat.NO.: PA01300

3th Edition

Description:The prion protein (Prp) is a protein of unknown function that is expressed in both normal cells and in cells affected by Transmissible spongiform encephalopathies (TSEs) which are lethal neurodegenerative diseases affecting numerous mammals.

Antigen:Synthetic peptide for N-terminal of Human Prion Protein

Form:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Confirmed by western blotting.