

**Rabbit IgG**  
**Polyclonal Antibody**  
**Cat.NO.: PA01308**

---

3th Edition

**Description:** The product is normal rabbit IgG purified from normal rabbit serum.

**Antigen:** The product is normal rabbit IgG purified from normal rabbit serum.

**Form:** Lyophilized product in PBS (without additives).

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**