

**mLST8 (L17) Anti-Human Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01310**

---

3th Edition

**Description:** mLST8 (mammalian lethal with sec-13), an about 39 kDa protein which has WD40 domains, has been identified as a configuration factor of mTORC1/mTOR complex 1 and then found to be also a configuration factor of mTORC2.

**Antigen:** Synthetic peptide of N-terminal of Human mLST8 (TAGYDHTVRFWQAHSGIC)

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**