

Anti-Mouse IgG (H+L) Goat IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01312

3th Edition

Description: Anti-Mouse IgG (H+L) Goat IgG Affinity Purify, 1 mg. Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃, purified with antigen.

Antigen: Mouse IgG (H+L)

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: The antibody is absorbed by rat IgG, goat IgG, rabbit IgG and human IgG. Cross-reactivity to goat IgG, rabbit IgG and human IgG is less than 0.1 %. Cross-reactivity to rat IgG is about 0.5 %.