

Anti-Mouse IgG (H+L) Rabbit IgG Fab' HRP

Polyclonal Antibody

Cat.NO.: PA01313

3th Edition

Description: Anti-Mouse IgG (H+L) Rabbit IgG Fab' HRP, 0.2 mg. Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃, purified with antigen peptide.

Antigen: Mouse IgG (H+L)

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: The antibody is absorbed by goat IgG, rabbit IgG and human IgG, but is not absorbed by rat IgG.