

Anti-Human/Mouse Phospho-FLT3 (Tyr842) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01698

3th Edition

Description: Receptor for the FL cytokine. Has a tyrosine-protein kinase activity.

Antigen: Synthesized peptide derived from human Flt-3 around the phosphorylation site of Y842.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Bone marrow cells.