

Anti-Human/Mouse/Rat Phospho-RAF1 (Ser301) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01702

3th Edition

Description:Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. Part of the Ras-dependent signaling pathway from receptors to the nucleus. Protects cells from apoptosis mediated by STK3.

Antigen:Synthesized peptide derived from human Raf-1 around the phosphorylation site of Ser301

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:In skeletal muscle, isoform 1 is more abundant than isoform 2.