

**Anti-Human/Mouse/Rat Phospho-DNM1 (Ser778) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01706**

---

3th Edition

**Description:** Microtubule-associated force-producing protein involved in producing microtubule bundles and able to bind and hydrolyze GTP. Most probably involved in vesicular trafficking processes. Involved in receptor-mediated endocytosis.

**Antigen:** Synthesized peptide derived from human Dynamin I around the phosphorylation site of Ser778

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**