

Anti-Human/Mouse/Rat Phospho-GRK2 (Ser685) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01709

3th Edition

Description: Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors, probably inducing a desensitization of them. Key regulator of LPAR1 signaling. Competes with RALA for binding to LPAR1 thus affecting the signaling properties of the receptor. Desensitizes LPAR1 and LPAR2 in a phosphorylation-independent manner.

Antigen: Synthesized peptide derived from human GRK 2 around the phosphorylation site of Ser685

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in peripheral blood leukocytes.