

**Anti-Human/Mouse Phospho-CREB2 (Ser219) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01719**

---

3th Edition

**Description:** Transcriptional activator. Binds the cAMP response element (CRE) (consensus: 5'-GTGACGT[AC][AG]-3'), a sequence present in many viral and cellular promoters. It binds to a Tax-responsive enhancer element in the long terminal repeat of HTLV-I.

**Antigen:** Synthesized peptide derived from human CREB-2 around the phosphorylation site of S219.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**