

**Anti-Human/Mouse/Rat Phospho-EPHA7 (Tyr791) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01722**

---

3th Edition

**Description:** Receptor for members of the ephrin-A family. Binds to ephrin-A1, -A2, -A3, -A4 and -A5.

**Antigen:** Synthesized peptide derived from human EphA7 around the phosphorylation site of Tyr791

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed.