

**Anti-Human/Mouse/Rat Phospho-FGFR1 (Tyr654) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01732**

---

3th Edition

**Description:** Receptor for basic fibroblast growth factor. Receptor for FGF23 in the presence of KL (By similarity). A shorter form of the receptor could be a receptor for FGF1 (aFGF).

**Antigen:** Synthesized peptide derived from human Flg around the phosphorylation site of Tyr654

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Detected in astrocytoma, neuroblastoma and adrenal cortex cell lines. Some isoforms are detected in foreskin fibroblast cell lines, however isoform 17, isoform 18 and isoform 19 are not detected in these cells.