

Anti-Human/Mouse/Rat Phospho-CD5 (Tyr453) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01737

3th Edition

Description: May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.

Antigen: Synthesized peptide derived from human CD5 around the phosphorylation site of Y453.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: