

Anti-Human/Mouse/Rat/Monkey Phospho-TH (Ser62) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01758

3th Edition

Description: Plays an important role in the physiology of adrenergic neurons.

Antigen: Synthesized peptide derived from human TH around the phosphorylation site of S62.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Mainly expressed in the brain and adrenal glands.