

**Anti-Human/Mouse/Rat/Monkey Phospho-TH (Ser62) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01758**

---

3th Edition

**Description:** Plays an important role in the physiology of adrenergic neurons.

**Antigen:** Synthesized peptide derived from human TH around the phosphorylation site of S62.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Mainly expressed in the brain and adrenal glands.