

Anti-Human/Mouse/Rat Phospho-CDK2 (Thr160) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01762

3th Edition

Description: Serine/threonine-protein kinase involved in the control of the cell cycle; essential for meiosis, but dispensable for mitosis.

Antigen: Synthesized peptide derived from human Cdk2/Cdc2 around the phosphorylation site of Thr160

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: