

Anti-Human/Mouse Phospho-KSR1 (Ser392) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01767

3th Edition

Description: Location-regulated scaffolding protein connecting MEK to RAF. Promotes MEK and RAF phosphorylation and activity through assembly of an activated signaling complex. By itself, it has no demonstrated kinase activity.

Antigen: Synthesized peptide derived from human Ksr-1 around the phosphorylation site of Ser392

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: