

Anti-Human/Mouse/Rat Phospho-DNM1 (Ser774) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01779

3th Edition

Description: Microtubule-associated force-producing protein involved in producing microtubule bundles and able to bind and hydrolyze GTP. Most probably involved in vesicular trafficking processes. Involved in receptor-mediated endocytosis.

Antigen: Synthesized peptide derived from human Dynamin I around the phosphorylation site of Ser774

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: