

Anti-Human/Mouse/Rat Phospho-BCL2 (Ser87) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01788

3th Edition

Description: Suppresses apoptosis in a variety of cell systems including factor-dependent lymphohematopoietic and neural cells. Regulates cell death by controlling the mitochondrial membrane permeability. Appears to function in a feedback loop system with caspases. Inhibits caspase activity either by preventing the release of cytochrome c from the mitochondria and/or by binding to the apoptosis-activating factor (APAF-1).

Antigen: Synthesized peptide derived from human Bcl-2 around the phosphorylation site of Ser87

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in a variety of tissues.