

Anti-Human/Mouse/Rat Phospho-P23 (Ser113) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01793

3th Edition

Description: Molecular chaperone that localizes to genomic response elements in a hormone-dependent manner and disrupts receptor-mediated transcriptional activation, by promoting disassembly of transcriptional regulatory complexes.

Antigen: Synthesized peptide derived from human p23 around the phosphorylation site of Ser113

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: