

**Anti-Mouse/Rat Phospho-TNNI3 (Ser22/S23) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01794**

---

3th Edition

**Description:** Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

**Antigen:** Synthesized peptide derived from human Troponin I-C around the phosphorylation site of Ser22/S23

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**