

Anti-Mouse/Rat Phospho-TNNI3 (Ser22/S23) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01794

3th Edition

Description: Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Antigen: Synthesized peptide derived from human Troponin I-C around the phosphorylation site of Ser22/S23

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: