

**Anti-Human/Mouse/Rat Phospho-XIAP (Ser87) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01798**

---

3th Edition

**Description:** Apoptotic suppressor. Has E3 ubiquitin-protein ligase activity. Mediates the proteasomal degradation of target proteins, such as caspase-3, SMAC or AIFM1. Inhibitor of caspase-3, -7 and -9. Mediates activation of MAP3K7/TAK1, leading to the activation of NF-kappa-B.

**Antigen:** Synthesized peptide derived from human XIAP around the phosphorylation site of Ser87

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous, except peripheral blood leukocytes.