

Anti-Human/Mouse/Rat Phospho-XIAP (Ser87) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01798

3th Edition

Description: Apoptotic suppressor. Has E3 ubiquitin-protein ligase activity. Mediates the proteasomal degradation of target proteins, such as caspase-3, SMAC or AIFM1. Inhibitor of caspase-3, -7 and -9. Mediates activation of MAP3K7/TAK1, leading to the activation of NF-kappa-B.

Antigen: Synthesized peptide derived from human XIAP around the phosphorylation site of Ser87

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous, except peripheral blood leukocytes.