

Anti-Human/Mouse/Rat Phospho-BRAF (Ser602) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01809

3th Edition

Description:Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron.

Antigen:Synthesized peptide derived from human Raf-B around the phosphorylation site of S602.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Brain and testis.