

**Anti-Human/Mouse/Rat Phospho-BRAF (Ser602) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01809**

---

3th Edition

**Description:**Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron.

**Antigen:**Synthesized peptide derived from human Raf-B around the phosphorylation site of S602.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Brain and testis.