

Anti-Human/Mouse Phospho-CD71 (Ser24) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01813

3th Edition

Description:CD71 is a cell surface receptor necessary for cellular iron uptake by the process of receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Mice that are deficient in this receptor show impaired erythroid development and abnormal iron homeostasis.

Antigen:Synthesized peptide derived from human CD71 around the phosphorylation site of S24.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: