

Anti-Human/Mouse/Rat Phospho-TPH1 (Ser58) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01818

3th Edition

Description: This gene encodes a member of the aromatic amino acid hydroxylase family. The encoded protein catalyzes the first and rate limiting step in the biosynthesis of serotonin, an important hormone and neurotransmitter. Mutations in this gene have been associated with an elevated risk for a variety of diseases and disorders, including schizophrenia, somatic anxiety, anger-related traits, bipolar disorder, suicidal behavior, addictions, and others.

Antigen: Synthesized peptide derived from human TPH1 around the phosphorylation site of S58.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Isoform 2 seems to be less widely expressed than isoform 1.