

**Anti-Human/Mouse/Rat/Monkey Phospho-CTTN (Tyr466) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01849**

---

3th Edition

**Description:**Contributes to the organization of the actin cytoskeleton and cell structure. Plays a role in the regulation of cell migration. Plays a role in the invasiveness of cancer cells, and the formation of metastases.

**Antigen:**Synthesized peptide derived from human Cortactin around the phosphorylation site of Tyr466

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**