

Anti-Human/Mouse/Rat/Monkey Phospho-CTTN (Tyr466) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01849

3th Edition

Description:Contributes to the organization of the actin cytoskeleton and cell structure. Plays a role in the regulation of cell migration. Plays a role in the invasiveness of cancer cells, and the formation of metastases.

Antigen:Synthesized peptide derived from human Cortactin around the phosphorylation site of Tyr466

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: