

Anti-Human/Mouse/Rat Phospho-c-Myc (Ser62) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01854

3th Edition

Description: Participates in the regulation of gene transcription. Binds DNA in a non-specific manner, yet also specifically recognizes the core sequence 5'-CAC[GA]TG-3'. Seems to activate the transcription of growth-related genes.

Antigen: Synthesized peptide derived from human c-Myc around the phosphorylation site of Ser62

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: