

Anti-Human/Mouse Phospho-ESPL1 (Ser801) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01860

3th Edition

Description: Caspase-like protease, which plays a central role in the chromosome segregation by cleaving the SCC1/RAD21 subunit of the cohesin complex at the onset of anaphase. During most of the cell cycle, it is inactivated by different mechanisms.

Antigen: Synthesized peptide derived from human Separase around the phosphorylation site of Ser801

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: