

Anti-Human/Mouse/Rat Phospho-GAB1 (Tyr659) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01865

3th Edition

Description:The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

Antigen:Synthesized peptide derived from human Gab 1 around the phosphorylation site of Y659.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: