

**Anti-Human/Mouse/Rat Phospho-p120 (Tyr228) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01870**

---

3th Edition

**Description:** Binds to and inhibits the transcriptional repressor ZBTB33, which may lead to activation of target genes of the Wnt signaling pathway (By similarity). May associate with and regulate the cell adhesion properties of both C- and E-cadherins. Implicated both in cell transformation by SRC and in ligand-induced receptor signaling through the EGF, PDGF, CSF-1 and ERBB2 receptors. Promotes GLIS2 C-terminal cleavage.

**Antigen:** Synthesized peptide derived from human p120 around the phosphorylation site of Tyr228

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in vascular endothelium.