

Anti-Human/Mouse/Rat Phospho-ERBB3 (Tyr1289) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01883

3th Edition

Description: Binds and is activated by neuregulins and NTAK.

Antigen: Synthesized peptide derived from human ErbB-3 around the phosphorylation site of Tyr1289

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Epithelial tissues and brain.