

Anti-Rat ADIPOQ Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01924

3th Edition

Description:Rat gAcrp30 (amp-1) is a 16.6 kDa protein that corresponds to the C-terminal globular domain of adiponectin (amp-1). gAcrp30/Adipolean is expressed and secreted exclusively by adipose tissue. This protein possesses unique pharmacological properties for regulation of body weight and lipid metabolism. It is a key component in controlling energy homeostasis. Recently, gAcrp30/Adipolean has been shown to inhibit the proliferation of myelomonocytic progenitors.

Antigen:Recombinant Rat Adiponectin protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Synthesized exclusively by adipocytes and secreted into plasma.