

**Anti-Mouse/Rat FABP2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01943**

---

3th Edition

**Description:**FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters. FABP2 is probably involved in triglyceride-rich lipoprotein synthesis. Binds saturated long-chain fatty acids with a high affinity, but binds with a lower affinity to unsaturated long-chain fatty acids. FABP2 may also help maintain energy homeostasis by functioning as a lipid sensor.

**Antigen:**Recombinant Mouse Fatty acid-binding protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in the small intestine and at much lower levels in the large intestine. Highest expression levels in the jejunum.