

Anti-Human/Rat CD48 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01953

3th Edition

Description: Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.

Antigen: Recombinant human CD48 antigen protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: