

Anti-Rat CD244 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01969

3th Edition

Description: Modulate other receptor-ligand interactions to enhance leukocyte activation.

Antigen: Recombinant human Natural killer cell receptor 2B4 protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in spleen, PBL, followed by lung, liver, testis and small intestine. Expressed not only in NK cells, but also on monocytes and basophils.