

**Anti-Human/Mouse UQCRFS1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01986**

---

3th Edition

**Description:**Component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is a respiratory chain that generates an electrochemical potential coupled to ATP synthesis.

**Antigen:**Recombinant human Cytochrome b-c1 complex subunit Rieske, mitochondrial protein

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**