

Anti-Human/Mouse UQCRFS1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01986

3th Edition

Description:Component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is a respiratory chain that generates an electrochemical potential coupled to ATP synthesis.

Antigen:Recombinant human Cytochrome b-c1 complex subunit Rieske, mitochondrial protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: