

Anti-Human/Mouse/Rat COX5B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01987

3th Edition

Description: This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

Antigen: Recombinant human Cytochrome c oxidase subunit 5B, mitochondrial protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: