

Anti-Human/Mouse TGM2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01991

3th Edition

Description: Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

Antigen: Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: