

Anti-Mouse BGLAP Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02005

3th Edition

Description: Constitutes 1-2% of the total bone protein. It binds strongly to apatite and calcium.

Antigen: Recombinant human BGLAP protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: