

**Anti-Human/Mouse/Rat FH Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02008**

---

3th Edition

**Description:** Fumarate hydratase, Also acts as a tumor suppressor.

**Antigen:** Recombinant human Fumarate hydratase protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in red blood cells; underexpressed in red blood cells (cytoplasm) of patients with hereditary non-spherocytic hemolytic anemia of unknown etiology.